

PART ONE: KEY CONCEPTS ON LASER THERAPY

- I. What is a laser?
- II. Laser therapy: **Interaction with the biological tissues**
- III. Laser therapy: **The therapeutic window**
- IV. Laser therapy: **The importance of average power**
- V. Laser therapy: **Photobiomodulation**
- VI. Laser therapy: **The effects of therapeutic laser light**

PART TWO: HOW TO MAKE PBM MORE EFFICIENT?

- I. PBM: **Which wavelengths can make it more efficient?**
- II. PBM: **Which emission mode can make it more efficient?**
- III. PBM: **Is a temperature control system necessary?**
- IV. PBM: **The correct dose**
- V. Summary: **How to make PBM more efficient?**

PART THREE : APPLICATION FIELDS

- I. Dental
- II. Podiatry
- III. Neurological
- IV. Dermatology
- V. Sport
- VI. Rehabilitation
- VII. Disc Pathologies
- VIII. Wounds
- IX. Wellness
- X. Acute trauma

PART FOUR : TREATMENT EXAMPLES

- I. Inflammations
- II. Nerve repair
- III. Œdema
- IV. Muscular injuries
- V. Tendon injuries
- VI. Pain management
- VII. Active and latent Trigger Points
- VIII. Spine pathologies

- IX. Treatment examples – regenerative medicine**
- X. Treatment examples – podiatry**
- XI. Treatment examples - wellness**
- XII. Treatment examples - otorhinolaryngology**

PART FIVE: MELECTRONIC LASER THERAPY DEVICES

- I. Laser line 2023
- II. MHEL
- III. THEAL

PART ONE: KEY CONCEPTS ON LASER THERAPY

I. What is **laser**?

Light

Amplification by the

Stimulated

Emission of

Radiation

Laser therapy is a medical treatment that uses laser light to stimulate a process called photobiomodulation (PBM)

I. What is **laser**?

PHYSICAL CHARACTERISTICS

Laser is an artificial electromagnetic radiation similar to the sun light, but laser light has some special features (monocromaticity, coherence, directionality) that allows to use it in some different applications.

I. What is **laser**?

LASER POINTERS



I. What is **laser**?

LASER METERS




I. What is **laser**?

SPEED CAMERAS



I. What is **laser**?

METAL LASER CUTTERS





I. What is **laser**?

LASER MARKING AND ENGRAVING



I. What is **laser**?

OPHTHALMOLOGY



I. What is **laser**?

SURGERY




I. What is **laser**?

DERMATOLOGY




I. What is **laser**?

DENTAL



I. What is **laser**?

REHABILITATION




I. What is **laser**?

Over 40 years of research and thousands of scientific researches have demonstrated the efficacy of laser light. Laser therapy is the physical therapy with **the most scientific publications**.


THERAPIES	SCIENTIFIC RESEARCH
LASER THERAPY	→ 91.225 ←
ULTRASOUND THERAPY	67.636
TENS	24.183
MICROWAVE THERAPY	7.377
SHOCKWAVE THERAPY	6.329
MAGNETOTHERAPY	4.245
TECAR	97

II. Laser therapy: Interaction with the biological tissues


What happens when laser light interacts with **biological tissues**?




II. Laser therapy: Interaction with the biological tissues




Reflection



Absorption




Transmission



Scattering

II. Laser therapy: Interaction with the biological tissues


A percentage of the laser light that hits the biological tissues is reflected.
Phototypes I and II reflect laser light more than phototypes V or VI.



II. Laser therapy: Interaction with the biological tissues

The percentage of reflected laser light must be reduced!

HOW?



II. Laser therapy: Interaction with the biological tissues


The percentage of reflected laser light must be reduced!

HOW?

- ❑ CLEANING THE PATIENT'S SKIN BEFORE TREATMENT
- ❑ PERFORMING THE TREATMENT ORTHOGONALLY TO THE TREATED TISSUE

II. Laser therapy: Interaction with the biological tissues


THE IMPORTANCE OF ORTHOGONALITY




II. Laser therapy: Interaction with the biological tissues

THE IMPORTANCE OF ORTHOGONALITY

The **scanning laser** allows to treat very large areas, but it is almost always **non-orthogonal** and the percentage of reflected laser light is higher.




II. Laser therapy: Interaction with the biological tissues



OUR FOCUS IS THE ABSORPTION COEFFICIENT THAT DETERMINES THE TRANSMISSION


II. Laser therapy: Interaction with the biological tissues

The laser light must interact with **the chromophores** in the human tissues, therefore the laser light that brings photon energy must be able to penetrate those tissues.



II. Laser therapy: Interaction with the biological tissues

The term **chromophore** defines a group of atoms capable of giving color to a substance.



II. Laser therapy: Interaction with the biological tissues

The laser light must penetrate inside the biological tissues. To do that, the laser wavelength must be not overly absorbed by the chromophores:




Melanin



Water



Hemoglobin



Hemoglobin oxide

III. Laser therapy: The therapeutic window

The laser light suitable for therapeutic activities/purposes must belong to the therapeutic window.

The therapeutic window is well defined in the scientific article:

“Mechanisms of low-level light therapy” (2006)

Proc. of SPIE Vol. 6140 614001-1

M. HAMBLIN, T. DEMIDOVA

III. Laser therapy: The therapeutic window

The therapeutic window is a range of wavelengths of laser light that can be used for therapeutic purposes.

“Mechanisms of low level light therapy” (2006)

Proc. of SPIE Vol. 6140 614001-1

M. HAMBLIN, T. DEMIDOVA


III. Laser therapy: The therapeutic window

“Mechanisms of low-level light therapy” (2006)


Proc. of SPIE Vol. 6140 614001-1

M. HAMBLIN, T. DEMIDOVA

**Therapeutic window
600nm - 1100nm**



III. Laser therapy: The therapeutic window



IV. Laser therapy: **The importance of average power**


The ability of a laser to penetrate inside biological tissues depends **solely on the wavelength**, but the number of photons and energy inside the tissues is proportional to the **average power applied: the greater the average power, the greater the energy carried in depth.**

IV. Laser therapy: **The importance of average power**

The ability of a laser to penetrate into biological tissues depends **solely on the wavelength**



IV. Laser therapy: The importance of average power



The higher the average power, the greater the depth of action of the laser beam

IV. Laser therapy: The importance of average power

The ability of a laser to penetrate into biological tissues depends **solely on the wavelength**



The peak power does not increase the ability of the laser light to penetrate inside biological tissues!

IV. Laser therapy: The importance of average power

Comparison of Light Penetration of Continuous Wave 810 nm and Superpulsed 904 nm Wavelength Light in Anesthetized Rats. Anders JJ¹, Wu X¹ - 2016 Sep;34(9):418-24. doi: 10.1089/pho.2016.4137. Epub 2016 Aug 8.

OBJECTIVE:

The purpose of this study was to investigate light transmission of continuous wave (CW) 810 nm wavelength light and 904 nm wavelength superpulsed light through skin and gastrocnemius muscle and skin only using an anesthetized Sprague-Dawley rat model.

RESULTS: The percentages of light transmission (fluence rate) through muscle and skin were

7.42% (810 nm wavelength)
4.01% (904 nm wavelength)

and through skin were

24.63% (810 nm wavelength)
19.94% (904 nm wavelength)

These data prove that transmission of CW 810 nm wavelength light through muscle and skin and skin alone is greater than transmission of superpulsed 904 nm wavelength light.

IV. Laser therapy: The importance of average power

High
Average power



High
Energy in depth



High
Depth of action



V. Laser therapy: **Photobiomodulation**

Laser therapy is a medical treatment that uses laser light to stimulate a process called **photobiomodulation (PBM)**


V. Laser therapy: **Photobiomodulation**

Photobiomodulation is defined as a form of light therapy that uses visible and near infrared laser sources.

It is a non-thermal or moderately thermal process that involves endogenous chromophores that cause photophysical and photochemical events at various biological scales. This process results in therapeutic outcomes not only aimed at reducing pain and inflammation, but also at promoting immunomodulation, wound healing and tissue regeneration.

V. Laser therapy: Photobiomodulation

There are over 90.000 scientific publications dedicated to laser therapy and over 7.000 dedicated to Photobiomodulation. In recent years, **Photobiomodulation is increasingly the subject of study and attention**, with an exponential increase in scientific publications produced annually. In fact, in the last 20 years we have gone from 80 researches published in 2000 to 1090 in 2020.



VI. Laser therapy: The effects of therapeutic laser light

The effects of therapeutic laser light on the tissue are three:

- **Photochemical**
- **Photothermal**
- **Photomechanical**

PHOTOCHEMICAL EFFECT

VI. Laser therapy: The effects of therapeutic laser light

PHOTOCHEMICAL EFFECT: LIGHT IS AT THE BASE OF LIFE

The light is essential for the wellness of living beings.



PLANTS



ANIMALS



HUMANS

VI. Laser therapy: The effects of therapeutic laser light

PHOTOCHEMICAL EFFECT: THE LIGHT EFFECT ON PLANTS



The bright energy is absorbed by a photoreceptor that triggers a fundamental chemical reaction for the wellness and the growth of the plants.

VI. Laser therapy: The effects of therapeutic laser light

PHOTOCHEMICAL EFFECT: THE LIGHT EFFECT ON PLANTS



Energy absorbed by the chlorophyll (photoreceptor)



Glucose + Oxygen

Essential sugar for
the wellness of plants

VI. Laser therapy: The effects of therapeutic laser light

PHOTOCHEMICAL EFFECT: THE LIGHT EFFECT ON PLANTS



Chlorophylline
photosynthesis

VI. Laser therapy: The effects of therapeutic laser light

PHOTOCHEMICAL EFFECT: THE LIGHT EFFECT ON HUMANS

There is a similar reaction in human beings.

The laser energy interacts with the human body to obtain a healthy and therapeutic effect.



Theal
THERAPY

VI. Laser therapy: The effects of therapeutic laser light

PHOTOCHEMICAL EFFECT: THE LIGHT EFFECT ON HUMANS

First law of photochemistry:

**“Light must be absorbed before
photochemistry can occur”**

VI. Laser therapy: The effects of therapeutic laser light

PHOTOCHEMICAL EFFECT: THE LIGHT EFFECT ON HUMANS


The photobiomodulation, at the base of metabolic activities, is triggered by the **selective and naturally thermal laser light absorption**, which provides three fundamental effects: **tissue regeneration, inflammatory action modulation and analgesic effect**.




VI. Laser therapy: The effects of therapeutic laser light

PHOTOCHEMICAL EFFECT:


The main photo acceptor that absorbs the bright energy to trigger photobiomodulation is the **Cytochrome C Oxidase** present in a cellular organelle: the mitochondrion.



CELL




MITOCHONDRION




CYTOCHROME

VI. Laser therapy: The effects of therapeutic laser light




When laser light is absorbed by cytochrome c it stimulates the electron transport chain to increase the production of adenosine triphosphate (ATP) within the mitochondria. Adenosine triphosphate (ATP) is the molecule that facilitates the transfer of energy within the cell.

VI. Laser therapy: The effects of therapeutic laser light




When the tissue is damaged, ATP production in the cell is impaired, which slows down the cell's metabolism as a protective mechanism. Photobiomodulation (PBM) helps restore the oxidative process which helps restore normal cell function.

VI. Laser therapy: The effects of therapeutic laser light




Laser stimulation also produces nitric oxide (NO). Nitric oxide (NO) is a powerful vasodilator and an important cell signaling molecule involved in many physiological processes. The release of nitric oxide (NO) increases circulation, decreases inflammation and improves the transport of oxygen and immune cells throughout the tissue.

VI. Laser therapy: The effects of therapeutic laser light



Laser stimulation also produces reactive oxygen species (ROS). Reactive oxygen species (ROS) have been shown to influence many important physiological signaling pathways, including the inflammatory response. Furthermore, the production of these signaling molecules has been shown to induce the **production of growth factors, increase cell proliferation and mobility and promote extracellular matrix deposition.**

VI. Laser therapy: The effects of therapeutic laser light



In concert, the production of these signaling molecules has been shown to:



Accelerate tissue regeneration


Modulate the inflammatory response

Reduce the pain sensation

PHOTOMECHANICAL EFFECT

VI. Laser therapy: The effects of therapeutic laser light

Laser light interacts with the tissue transforming the light energy in an acoustic / mechanical impulse.



VI. Laser therapy: The effects of therapeutic laser light

Interaction with mechanoreceptors



Analgesic effect



Mechanical stress of tissues





Supports tissue regeneration

Anti-inflammatory action



VI. Laser therapy: The effects of therapeutic laser light


980 nm: Mainly absorbed by water, it allows to optimize the action on thermoreceptors and mechanoreceptors.



VI. Laser therapy: The effects of therapeutic laser light


980 nm

When combined with the E²C mode, ensure the right interaction with the peripheral nervous system, activating the Gate Control mechanism for a quick analgesic effect.




VI. Laser therapy: The effects of therapeutic laser light

1064 nm: Its high directionality allows to convey the correct dose of energy directly to the noxa.



"Mechanisms of low level light therapy" (2006)
Proc. of SPIE Vol. 6140 614001-1 - M. HAMBLIN, T. DEMIDOVA




"Re-evaluation of the near infrared spectra of mitochondrial cytochrome c oxidase: Implications for non invasive in vivo monitoring of tissues"
(2014) BBA Bioenergetics - G.MASON, P. NICHOLLS, E. COOPER

VI. Laser therapy: The effects of therapeutic laser light

1064 nm

The result is a perfect synergy that harmonizes the rapid analgesic effect with a control of inflammatory processes and the deep activation of vital metabolic processes for all cellular activities.



VI. Laser therapy: The effects of therapeutic laser light

980 nm and 1064 nm

Due to these peculiarities, these two wavelengths trigger further metabolic pathways that may act in some cases at the same time as those triggered by photobiomodulation.

VI. Laser therapy: The effects of therapeutic laser light



PHOTOTHERMAL EFFECT


VI. Laser therapy: The effects of therapeutic laser light



The laser light interacts with the tissue transforming the light energy in heat.

VI. Laser therapy: The effects of therapeutic laser light


The thermal gradient created increases the blood flow at the local level.



VI. Laser therapy: The effects of therapeutic laser light

All the wavelengths of the therapeutic window applied with a sufficiently high-power density allow to obtain a photothermal effect.

Among these the 980nm, having a peak of absorption on the water, at the same power can create a greater thermal gradient.




"Mechanisms of low-level light therapy" (2006)
Proc. of SPIE Vol. 6140 614001-1 - M. HAMBLIN, T. DEMIDOVA


VI. Laser therapy: The effects of therapeutic laser light

Effects associated with the photothermal action:

Decontracting




Anti-edema



Muscle relaxant



PART TWO: HOW TO MAKE PBM MORE EFFICIENT?



● Wavelengths

● Thermal control

● Emission mode

● High Average Power

● Correct transfer method

I. PBM: Which wavelengths can make it more efficient?

Does the **Cytochrome C Oxidase** absorb all the wavelengths of the therapeutic window in the same way?

NO!



I. PBM: Which wavelengths can make it more efficient?


“Multiple Roles of Cytochrome c Oxidase in Mammalian Cells Under Action of Red and IR-A Radiation» (2010) IUBMB LIFE - T. I. KARU



There are numerous studies that explain **Photobiomodulation**. **Tina Karu** is one of the most important expert and researcher of the phenomenon of light absorption on **Cytochrome C Oxidase**.

I. PBM: Which wavelengths can make it more efficient?

Several studies show how **Cytochrome C Oxidase** mainly absorbs laser light:




"Re-evaluation of the near infrared spectra of mitochondrial cytochrome c oxidase: Implications for non invasive in vivo monitoring of tissues"

BBA Bioenergetics - (2014) G.MASON, P. NICHOLLS, E. COOPER

I. PBM: Which wavelengths can make it more efficient?

Numerous scientific articles have shown how the wavelengths in the range **600 nm - 670 nm** and in the range **780 nm - 880 nm** are more absorbed by the main **photoacceptor** of laser therapy (Cytochrome C Oxidase) and therefore make it possible to make **photobiomodulation** (PBM) more efficient.

I. PBM: Which wavelengths can make it more efficient?



I. PBM: Which wavelengths can make it more efficient?

Other wavelengths within the therapeutic window (for example 980 nm and 1064 nm), although less absorbed by the Cytochrome C oxidase, are useful because they guarantee excellent interaction with thermo- and mechanoreceptors.

- **980 nm: the most absorbed by water within the therapeutic window**
- **1064 nm: less scattering and more directional**

II. PBM: Which emission mode can make it more efficient?

It has been demonstrated that the optimal laser emission to trigger photobiomodulation must be continuous or pulsed with a pulse duration of at least a few milliseconds (10^{-3} seconds), as stated by Hamblin, in 2010 in the article “Effect of Pulsing in Low-Level Light Therapy”.


II. PBM: Which emission mode can make it more efficient?



ATTENTION !!!
SUPERPULSED LASER: pulse duration microseconds ($s\ 10^{-6}$)
nanoseconds ($s\ 10^{-9}$)



Pulse duration



II. PBM: Which emission mode can make it more efficient?

CLINICAL TEST

Key words: Hifitrap®; pulsed Nd:YAG laser; photomechanical effect; tissue repair; extracellular matrix

Living for Health [05]

Relationship between cellular and systemic effects of pulsed Nd:YAG laser.

Gialdai F¹., Monici M².,

1) AS/Ampur Joint Laboratory, Dept. Clinical Physiopathology, Florence University, Florence, Italy;
2) AS/Ampur Joint Laboratory, ASA Bio. Div., Dept. Clinical Physiopathology, Florence University, Florence, Italy.

ABSTRACT

Notwithstanding the wide diffusion of lasertherapy in clinics and numerous studies reported in literature, molecular mechanisms of interaction between laser and tissues are not well understood.

The analysis of biological effects induced by laser radiation is rather complicated due to the wide possibility of setting instruments, the variability of applied protocols and the differences in treated tissues.

In this review, we describe our studies on the cellular and molecular mechanisms at the basis of the systemic effects produced by treatment with pulsed Nd:YAG laser, that is known as Hifitrap.

Starting from studies on photothermal effects, the hypothesis is that this type of laser cause an indirect photomechanical effect. The heat produced by transfer of radiation energy to the irradiated volume, diffuse into surrounding tissues, inducing temperature gradients which result in transitory modifications of mechanical-elastic properties of the extracellular microenvironment. These changing mechanical forces acting on cells.

Considering these studies and knowing

the key role of the extracellular matrix, not only as a structural support but also in maintaining tissue homeostasis, our experiments focused on the analysis of extracellular matrix molecules and cytoskeleton behavior, responsible of contact between cell and matrix and considered the best candidate to act as a mechanotransducer.

The data obtained have shown, in laser-treated cells, an increase in production of ECM molecules, such as aggrecan, collagen I and II, and a reorganization of microtubules and actin microfilaments network. It is well known that similar effects are obtained when cells are subjected to mechanical stress. Our data on absorption of Nd:YAG pulses by matrix components (proteins and polysaccharides) suggest that Nd:YAG pulses principally interact with the extracellular matrix, whose transitory deformation applies a mechanical stress to the cells.

We then focused on the effect of pulsed Nd:YAG on endothelial function and tissue repair processes. In treated endothelial cells and fibroblasts, key elements of angiogenesis and tissue repair, we found

overexpression of genes involved in the chemokine-mediated inflammatory pathways. Moreover, the treatment promoted the formation of ordered endothelial monolayers as well as ordered fibronectin fibril assembling. The findings indicate that treatment with Nd:YAG pulses has a stimulatory effect in the acute phase of inflammation and significant effect on the remodeling phase of tissue repair, also considering the important role that fibronectin plays in tissue structure regeneration. Therefore we can support that Hifitrap can efficaciously promote tissue repair processes.

INTRODUCTION

In spite of a wide application in clinics, many studies and a great body of literature, the molecular mechanisms of the interaction between laser and tissues, and the consequent cellular response, are still not completely known. They are object of current and future research in the field of laser biomedical application.

Unfortunately, not always scientifically rigorous studies, a limited knowledge of the molecular and cellular mechanisms underlying the biological effects of laser and, in turn, the systemic effects of laser therapy give rise to controlling results, unsupported hypotheses and unconvincing theories.

The studies on laser biological effects are very difficult due to the variety of biological responses that depend on laser source (wavelength, continuous/pulsed mode), operative conditions (fluence, time of exposure, etc.) and biological substrate considered (the body area, the tissue, the cell type etc.). Nevertheless, they are of critical importance for correct clinical applications, to improve instruments and protocols, to increase therapeutic effectiveness.

When the light interacts with a biological tissue a small part of radiation (~3-5%) is specularly reflected, the most part propagates within the tissue and it is partially diffused (scattering) and

“Relationship between cellular and systemic effects of pulsed Nd:YAG laser.”

CIALDAI F., MONICI M.

Starting from studies on photothermal effects, the hypothesis is that this type of laser cause an indirect photomechanical effect. The heat produced by transfer of radiation energy to the irradiated volume, diffuse into surrounding tissues, inducing temperature gradients which result in transitory modifications of mechanical-elastic properties of the extracellular microenvironment, thus changing mechanical forces acting on cells.

II. PBM: Which emission mode can make it more efficient?

CLINICAL TEST

Key words: Photomechanical stress, pulsed Nd:YAG laser, connective tissue, extracellular matrix, Hilterapia®

Effects of pulsed Nd: YAG laser at molecular and cellular level. A study on the basis of Hilterapia®.

Monici M., Cialdai F., Fusi F., Romano G., Pratesi R.¹

1. Whetstone MS Research Division, Dept of Clinical Pharmacology, University of Florence, Italy

2. Dept. of Clinical Pharmacology, University of Florence, Via Pieraccini 6, 50139 Florence, Italy

3. Dept. of Physics, University of Florence, Via Sansone 1, I-50197 Santa Lucia, Florence, Italy

ABSTRACT

Lasers have been widely applied in many different fields of medicine, proving their effectiveness in the treatment of a wide range of diseases. In spite of the great amount of literature, it is difficult to understand the molecular and cellular mechanisms at the basis of the systemic effects induced by laser irradiation because of different kinds of laser used, operative conditions, variety of biological targets and responses.

The application of high power lasers in physiotherapy is quite recent. It is due to the development of instruments which allow the control of photothermal and photochemical processes, so as to obtain therapeutic effects without tissue damage. In particular, pulsed Nd:YAG laser has proved feasibility and efficacy in the treatment of many different musculoskeletal diseases and it is believed to have anti-inflammatory, antioedema, analgesic and also reparative effects.

The aim of the studies here presented was to contribute in understanding the molecular mechanisms and cellular processes at the basis of the systemic effects produced by pulsed Nd:YAG laser irradiation.

Owing to the lack of chromophores efficiently absorbing Nd:YAG radiation

(wavelength 1064 nm) in cells and tissues, we hypothesized that, rather than photochemical processes, aspecific mechanisms probably due to combined photothermal and photomechanical interactions could be responsible for the above mentioned effects of pulsed Nd:YAG laser.

The finding suggests that cells "sense" pulsed Nd:YAG laser irradiation and respond to it through mechanotransduction machinery. We hypothesize that the interaction between tissue and laser radiation alters the mechanism of cell microenvironment, thus acting on the cells as a mechanical stress.

INTRODUCTION

Phototherapy, that is the use of light for the treatment and prevention of diseases, has been widely used from ancient times till now. From the time of the Pharaohs until relatively recent times the source of light was the sun.

The last century saw a rapid evolution in light sources, from inefficient arc lamps to lasers, which are the most advanced kind of light source.

The great advantage of the laser, in comparison with other sources, is the very high intensity and monochromaticity of the emitted radiation and also the

possibility to be effectively focused and coupled to optical fibres.

Lasers have been widely applied in many different fields of medicine, proving their effectiveness in the treatment of a wide range of diseases [1, 2].

In spite of the great amount of literature, the molecular and cellular mechanisms at the basis of the systemic effects induced by laser irradiation are mostly unknown. The studies on this subject are very difficult because of the numerous effects and the variety of biological responses, the kind of laser used, the operative conditions, the biological targets (different areas of the body, different tissues, different cell populations, etc. ...). However, they are very important because the increase in knowledge can lead to a higher therapeutic efficacy by improvement of laser sources and treatment protocols.

Depending on interaction time and effective power density, three types of interactions between laser radiation and tissues can be distinguished: photothermal, photochemical and photomechanical [3]. The effects induced by low power lasers, the first to be applied in physiotherapy, are mostly due to photochemical processes. These occur when endogenous or exogenous chromophores introduced in the tissue absorb radiation of suitable wavelength.

A chromophore molecule which absorbs a photon is converted in an excited state and may subsequently participate in a chemical reaction that leads to the final biological effect [4].

High power lasers have been used at first for tissue ablation and surgery, because they are able to produce important photothermal and photomechanical effects (stress waves) [5]. Their application in fields different from surgery, such as physiotherapy, is quite recent and it has been possible thanks to the development of laser systems with emission modalities which allow the control of photothermal and photomechanical processes, so as to obtain therapeutic effects without tissue damage.

In particular, pulsed Nd:YAG laser has proved its versatility and efficacy in the

“Effects of pulsed Nd: YAG laser at molecular and cellular level. A study on the basis of Hilterapia.”

MONICI M., CIALDAI F., FUSI F., ROMANO G., PRATESI R.


Owing to the lack of chromophores efficiently absorbing Nd:YAG radiation (wavelength 1064 nm) in cells and tissues, we hypothesized that, rather than photochemical processes, aspecific mechanisms probably due to combined photothermal and photomechanical interactions could be responsible for the above mentioned effects of pulsed Nd:YAG laser.

In conclusion, our results demonstrate that the effects of Nd:YAG pulses on culture cells are very similar to those induced by mechanical stress, thus supporting our

III. PBM: Is a temperature control system necessary?

Photobiomodulation is a **non-thermal** or **moderately thermal** process. It is essential to monitor the temperature of biological tissues during laser therapy and modulate the laser therapy according to the thermal response of biological tissues (**Vobulate Thermal Control**).

III. PBM: Is a temperature control system necessary?



III. PBM: Is a temperature control system necessary?



III. PBM: Is a temperature control system necessary?

- Treatment of pathologies in the acute phase
- Control of the excessive vascularization of the transition tissues to limit the effects of exacerbation or avoid possible blood effusions in case of micro lesions
- Treatment of neuropathies
- Greater safety in the treatment of dark skin types and tattoos

IV. PBM: The correct dose



PHOTOMEDICINE AND LASER
SURGERY

Mary Ann Liebert, Inc. publishers

Journals

Search

Alerts

[Photomed Laser Surg.](#) 2013 May; 31(5): 189–191.

doi: [10.1089/pho.2013.3510](#)

PMCID: PMC3643261

PMID: [23600376](#)

Is It Time to Consider Photobiomodulation As a Drug Equivalent?

[Tiina Karu](#), PhD, DrSci^M

► Author information ► Copyright and License information ► [Disclaimer](#)

This article has been [cited by](#) other articles in PMC.

THE QUESTION OF WHETHER photobiomodulation should be used as a drug equivalent arose in my mind after listening to presentations at the recent conference of the World Association for Laser Therapy (WALT)-2012 (Gold Coast City, Australia), and later at home when searching MEDLINE® for the years 2009–2012. Photobiomodulation (earlier terms: low level laser therapy, LLLT, laser biostimulation) has been used in clinical practice for >40 years by now, and its action mechanisms on cellular and molecular levels have been studied for >30 years. Enthusiastic medical specialists successfully used photobiomodulation in treating healing-resistant wounds and ulcers (e.g., chronic diabetic ulcers), in pain management, and in spinal cord and nervous system injuries when other methods had had limited success.¹ However, photobiomodulation is still not a part of mainstream medicine. The goal of the present Editorial is to highlight some important recent developments in clinical applications and in studies of cellular and molecular mechanisms behind the clinical findings.

Is it Time to Consider Photobiomodulation As a Drug Equivalent?

Photomedicine and Laser Surgery (2013)

T.Karu

IV. PBM: The correct dose

The **dose** is the determined quantity of a substance, in relation to an effect to be achieved directly or through the presence of other substances in reciprocal quantitative relationship. In the rehabilitation field, we mean **the amount of energy needed to induce cellular metabolic activity**.

The energy dose can be insufficient, effective or toxic.

IV. PBM: The correct dose

In order to precisely define the calculation of the therapeutic dose we must consider the following physical quantities:

- 1) **Average Power (W)**
- 2) **Laser Beam Area (cm²)**

$$(\text{Average Power} / \text{Laser Beam Area}) = \text{POWER DENSITY (W/ cm}^2\text{)}$$

IV. PBM: The correct dose




$$\text{POWER DENSITY (IRRADIANCE)} = \frac{\text{AVERAGE POWER (W)}}{\text{LASER BEAM AREA (cm}^2\text{)}} = \frac{10 \text{ W}}{10 \text{ cm}^2} = 1 \text{ W/cm}^2$$

IV. PBM: The correct dose



$$\text{POWER DENSITY (IRRADIANCE)} = \frac{\text{AVERAGE POWER (W)}}{\text{LASER BEAM AREA (cm}^2\text{)}} = \frac{10 \text{ W}}{20 \text{ cm}^2} = \mathbf{0.5 \text{ W/cm}^2}$$

IV. PBM: The correct dose




DOSE = POWER DENSITY X EXPOSURE TIME

$$\frac{\text{AVERAGE POWER (W)}}{\text{LASER BEAM AREA (cm}^2\text{)}} \times \text{EXPOSURE TIME (sec)} = \frac{\text{J}}{\text{sec} \times \text{cm}^2} \times \text{sec} = \frac{\text{J}}{\text{cm}^2} = 10 \frac{\text{J}}{\text{cm}^2}$$

$$\text{APPLIED ENERGY} = \text{DOSE} \times \text{AREA} = \frac{\text{J}}{\text{cm}^2} \times \text{cm}^2 = 100 \text{ J}$$

IV. PBM: The correct dose



Photobiomodulation delivery parameters: an evidenced based approach

Photobiomodulation, photomedicine and Laser Surgery (2021)


M.Grootveld, M. Cronshaw and S. Parker

Feeling the Heat: Evolutionary and Microbial Basis for the Analgesic Mechanisms of Photobiomodulation Therapy

Photobiomodulation, photomedicine and Laser Surgery (2019)

M. Cronshaw, S. Parker and P. Arany

IV. PBM: The correct dose



HEALING

Superficial targets: 2 – 10 J/cm² [best 5 J/cm²]

Deeper structures (to 1cm): 20 – 100 J/cm²
[best 50 J/cm²]


ANALGESIA AND ANTI-INFLAMMATORY

Superficial targets: 10 – 30 J/cm²

Deeper structures (to 1cm): 100 – 300 J/cm²


V. Summary: How to make PBM more efficient?

1. CORRECT TRANSFER METHOD (Orthogonality)



V. Summary: How to make PBM more efficient?

2. **WAVELENGTHS** in the range **600 nm - 670 nm** and **780 nm - 880 nm**



V. Summary: How to make PBM more efficient?

2. WAVELENGTHS in the range **600 nm - 670 nm** and **780 nm - 880 nm**

Wavelengths in the range **900 nm - 1100 nm** have an excellent interaction with thermo and mechanoreceptors.

- **4 WAVELENGTHS: 650 nm + 810 nm + 980 nm + 1064 nm**
- **Possibility to activate and deactivate the wavelengths**
- **Ability to create custom mixes**

V. Summary: How to make PBM more efficient?

3. Emission mode: **Continuous** or **pulsed emission** with a pulse duration of at least a few milliseconds ($s\ 10^{-3}$)



ATTENTION !!!
SUPERPULSED LASER: pulse duration microseconds ($s\ 10^{-6}$)
nanoseconds ($s\ 10^{-9}$)



V. Summary: How to make PBM more efficient?


4. **THERMAL CONTROL:** it is not sufficient a thermal control, but it is necessary a precise and powerful thermal control. **Photobiomodulation** is a **non-thermal** or **moderately thermal** process. It is essential to monitor the temperature of biological tissues during laser therapy and modulate the laser therapy according to the thermal response of biological tissues (**Vobulate Thermal Control**).

V. Summary: How to make PBM more efficient?

5. HIGH AVERAGE POWER



High Average Power to bring the right dose to the biological tissue



● Wavelengths

Wavelengths in the range
600 nm - 670 nm and
780 nm - 880 nm

● Emission mode

Continuous or **pulsed emission** with a pulse duration of at least a few milliseconds

● Correct transfer method

Reduce energy losses to
always transfer the **correct dose** (orthogonal treatment)

● Thermal control

Vobulate thermal control

● High Average Power

PART THREE : APPLICATION FIELDS



DENTAL

- Herpes labialis
- Angular cheilitis
- Dental hygiene
- Ulcers
- Leukoplakia
- Mouth therapy
- ATM
- Conservative therapy
- Bactericide action



DENTAL

- Herpes labialis
- Angular cheilitis
- Dental hygiene
- Ulcers
- Leukoplakia
- **Mouth therapy**
- ATM
- Conservative therapy
- Bactericide action



PODIATRY

- Onychomycosis
- Plantar warts
- Diabetic ulcers
- Vascular malleolus ulcers
- Foot pathologies




PODIATRY

- Onychomycosis
- Plantar warts
- Diabetic ulcers
- Vascular malleolus ulcers
- **Foot pathologies**



NEUROLOGICAL



DERMATOLOGY



SPORT

- Pain therapy
- Traumatology
- Accident prevention
- Accelerating physical recovery
- Fast elimination of edemas
- Physiological accelerator
- Muscular lesions
- Acute pathologies
- Tendinopathy



SPORT

- Pain therapy
- Traumatology
- Accident prevention
- Accelerating physical recovery
- Fast elimination of edemas
- Physiological accelerator
- Muscular lesions
- Acute pathologies
- **Tendinopathy**



REHABILITATION

- Pain therapy
- Acute and chronic pathologies
- Increased mobility
- Lymphatic drainage
- Traumatology
- Rheumatology
- Tissue repair and regeneration



REHABILITATION

- Pain therapy
- Acute and chronic pathologies
- Increased mobility
- Lymphatic drainage
- Traumatology
- Rheumatology
- Tissue repair and regeneration



DISC PATHOLOGIES



DISC PATHOLOGIES



WOUNDS

- Wounds
- Diabetic ulcers
- Vascular ulcers
- Decubitus ulcers and lesions



WOUNDS

- Wounds
- Diabetic ulcers
- Vascular ulcers
- Decubitus ulcers and lesions



WELLNESS



ACUTE TRAUMA



ACUTE TRAUMA

PART FOUR: TREATMENTS EXAMPLES

Treatment examples

INFLAMMATIONS

NERVE REPAIR

ŒDEMAS

MUSCLE INJURIES

TENDON INJURIES

PAIN MANAGEMENT

ACTIVE AND LATENT TRIGGER POINTS


SPINE PATHOLOGIES

I. Inflammations

- PBM activates microcirculation
- PBM promotes the reduction of inflammation
- PBM breaks the cycle of chronic inflammation

Strengths of Mechatronic Medical lasers

1. The best range of wavelengths, emission modes and power to maximise PBM
2. Temperature control system to treat acute inflammation and avoid pain return (THEAL and CHELT)
3. Different applicators to adapt the therapy to the treatment area and to transfer the correct dose in orthogonal mode. In addition to the wide range of powers, it is possible to deliver the right dose at different depths depending on the depth of the pathology.




II. Nerve repair

Injuries such as cuts, diabetes, sprains, carpal tunnel syndrome

- PBM not only treats the symptoms, but repairs the nerves.
- Improves nerve function and reduces inflammation- PBM can avoid surgery or improve results after surgery.

Strengths of Mectronic Medical lasers

1. The best range of wavelengths, emission modes and power to maximise PBM
2. Temperature control system to treat acute inflammation and avoid pain return (THEAL and CHELT)
3. Different applicators to adapt the therapy to the treatment area and to transfer the correct dose in orthogonal mode. In addition to the wide range of powers, it is possible to deliver the right dose at different depths depending on the depth of the pathology.



II. Nerve repair

Treatment of carpal tunnel with THEAL Therapy



III. Oedema

- PBM promotes vasodilation and activates microcirculation by increasing NO production (also without heat)
- PBM reduces pain and oedema after injuries and surgeries.

Strengths of Mectronic Medicales lasers

1. Best wavelengths, emission modes and power
2. Temperature control system to act on the anti-edema effect in the acute phase (THEAL e CHELT)
3. Different applicators to adapt the therapy to the treatment area and to transfer the correct dose in orthogonal mode.
4. Wide range of power: it is possible to use the right dose at different depths depending on the pathology
5. Pre-defined protocols

Use of low-level laser therapy to reduce postoperative pain, edema, and trismus following third molar surgery: A systematic review and meta-analysis

Francisca Jemifer Duarte de Oliveira ¹, Giuliana Moura Luz Cordeiro Brasil ²,
Gabriella Peixoto Araujo Soares ³, Daniel Felipe Fernandes Paiva ⁴,
Francisco de Assis de Souza Junior ⁵

Affiliations: [+ expand](#)

PMID: 34217567 · DOI: 10.1016/j.joms.2021.06.006

Abstract

This Systematic Review and Meta-analysis was conducted with the following PICO question: patients undergoing third molar surgery (P) can benefit from low-intensity laser therapy (I) as compared to other postoperative management (C) to reduce pain, edema, and trismus (O), evaluated in previous randomized clinical trials (S). Databases used were PubMed, SCOPUS, Web of Science, and Biblioteca Virtual em Saúde, screening for studies published between 2015 and 2020. The meta-analysis was based on the standardized mean difference (SMD), under a 95% confidence interval (CI). 246 studies were initially included, and after the screening of data, 10 studies were selected for the final sample. The qualitative analysis resulted in favorable results for pain and edema management in most studies, whereas trismus remained controversial. Meta-analysis resulted in (SMD, -0.53; 95% CI, -0.82, -0.24), (SMD, -0.60; 95% CI, -0.81, -0.39), and (SMD, -0.62; 95% CI, -2.63, 1.39) for pain, edema, and trismus, respectively, indicating statistical success on pain and edema reduction, but not for trismus. LLT can act on reducing postoperative pain and edema following third molar surgery, whereas, trismus remains not significantly changed.

III. Oedema




BEFORE



AFTER

IV. Muscular injuries

- Common muscle injuries, especially in sports
- Typical pain on contraction, stretching and direct pressure
- Post-surgical muscle trauma
- PBM reduces inflammation
- PBM reduces muscle fatigue
- PBM reduces scarring
- PBM enables stronger repair and reduces the possibility of relapse



Strengths of Mectronic Medicales lasers

1. Best wavelengths, emission modes and power
2. Temperature control system to act on the anti-edema effect in the acute phase (THEAL e CHELT)
3. Different applicators to adapt the therapy to the treatment area and to transfer the correct dose in orthogonal mode.
4. Wide range of power: it is possible to use the right dose at different depths depending on the pathology
5. Pre-defined protocols

IV. Muscular injuries



IV. Muscular injuries



V. Tendon injuries

- Effective for tendonitis and also for tendinosis
- Tendonitis: acute and highly inflamed
- Tendinosis: a degenerative and chronic condition
- PBM improves inflammation and collagen (type 1) by stimulating cells
- PBM allows scar tissue to shrink, repair itself more firmly and be less susceptible to injury

Strengths of Mectronic Medicales lasers

1. The best wavelengths, emission modes and power
2. Temperature control system to treat acute inflammation (tendonitis) and chronic inflammation (tendinosis) (THEAL and CHELT)
3. Different applicators to adapt the therapy to the treatment area and to transfer the correct dose in orthogonal mode.
4. Wide power range to deliver the right dose at different depths depending on the pathology.
5. Predefined protocols



Achilles tendon treatment with THEAL Therapy



VI. Pain Management

- PBM activates microcirculation
- PBM promotes reduction of inflammation
- Photomechanical effect of PBM for peripheral pain information

Strengths of Mectronic Medicales lasers

1. The best wavelengths, emission modes and power
2. Temperature control system to activate the analgesic effect through heatless photomechanical laser effect (THEAL and CHELT)
3. Patented E2C stochastic emission mode for variable photomechanical effect to obtain pain information (analgesic effect)
4. Different applicators to adapt the therapy to the treatment area and deliver the correct dose in orthogonal mode
5. Pre-defined protocols

Effectiveness of High Intensity Laser Therapy for Reduction of Pain in Knee Osteoarthritis

Anna Angelova¹, Elena Mileva¹

Affiliations: + expand

PMID: 28096711 PMCID: PMC5206453 DOI: 10.1155/2016/9163618

[Free PMC article](#)

Abstract

Introduction. Osteoarthritis is the most common type of arthritis. It is the main cause of chronic musculoskeletal pain and disability among the elderly population. **Aim.** This is a pilot, randomized clinical study about the effect of high intensity laser therapy in patients with osteoarthritis of the knee (OA of the knee). **Material and Method.** 72 patients (aged between 39 and 83 years) with (clinically and radiographically proved) OA of the knee were included in the study. They were randomized in two groups: therapeutic (test) one (n = 37, 65.11 ± 1.40 (mean ± SD) years old; patients were treated with HILT) and control group (n = 35, 64.71 ± 1.98; patients receive sham laser). Both groups had seven sessions of treatment. VAS and dolorimetry were used for assessment of pain before and after the therapy. Pedobarometric analysis (static and dynamic) was used to assess comparatively the contact surface area and maximum pressure under the heel. **Results.** Pain levels measured by VAS and dolorimetry decreased significantly in the therapeutic group after seven days of treatment (p < 0.001). **Conclusion.** The results after seven days of treatment show more intensive and cumulative effect after the application of high intensity laser therapy in comparison to sham laser. This is the reason why HILT can be a method of choice in the treatment of gonarthrosis.

VII. Active and latent Trigger Points

- PBM activates microcirculation
- PBM promotes reduction of inflammation
- Photomechanical effect of PBM for peripheral pain information

Strengths of Mectronic Medicales lasers

1. The best wavelengths, emission modes and power
2. Temperature control system to activate the analgesic effect through heatless photomechanical laser effect (THEAL and CHELT)
3. Patented E2C stochastic emission mode for variable photomechanical effect to obtain pain information (analgesic effect)
4. Different applicators to adapt the therapy to the treatment area and deliver the correct dose in orthogonal mode
5. Pre-defined protocols

Comparison of the effects between lasers applied to myofascial trigger points and to classical acupoints for patients with cervical myofascial pain syndrome

Wen-Han Chang¹, Li-Wen Tu², Yu-Cheng Pui³, Chih-Kuang Chen³, Szu-Hung Wang⁴, Alice-Ms Wong⁵

Affiliations: + expand

PMID: 35166212 PMCID: PMC8847842 DOI: 10.1016/j.tj.2020.05.020

Free PMC article

Abstract

Background: To compare the immediate effectiveness of low-level laser therapy (LLLT) applied to classical acupoints versus trigger points for patients with cervical myofascial pain syndrome (MPS).

Methods: This was a single-blinded, randomized, placebo-controlled trial. This study was performed in a university-affiliated medical center. One hundred participants with cervical myofascial pain syndrome were randomly allocated to four treatment groups, including (1) acupoint therapy (AcuT), (2) acupoint control (AcuC), (3) trigger point therapy (TriT), and (4) trigger point control (TriC) groups. Low-level laser (810-nm) therapy was used in both therapy groups, while the same procedure was performed without laser in the acupoint control groups. The patients were evaluated based on visual analogue scale (VAS) pain score, pressure pain threshold, and cervical range of motion (ROM) before and after the therapy.

Results: Immediate pain relief was observed in the TriT group ($p < 0.01$). The TriT group showed improved cervical ROM in ipsilateral bending ($p < 0.01$), while the AcuT group did not.

Conclusions: LLLT applied to trigger points could significantly relieve myofascial pain and was effective in relieving cervical ROM limitations. Considering the risk of pneumothorax, laser therapy at trigger points for patients with cervical MPS may be a choice when acupuncture therapy is unavailable.

VII. Active and latent Trigger Points

Treatment of an active Trigger Point with THEAL Therapy



VII. Active and latent Trigger Points

Treatment of an active Trigger Point with THEAL Therapy




VIII. Spine pathologies

- Lumbago, lumbosciatica, cervico-brachialgia, stenosis, herniated disc, etc...
- PBM activates the microcirculation
- PBM promotes the modulation of inflammation
- PBM reduces compression and stenosis of the spinal canal

Strengths of Mectronic Medicales lasers

1. The best range of wavelengths, emission modes and power
2. Temperature control system to treat acute nerve inflammation with a special collimated applicator (THEAL e CHELT)
3. Patented E2C stochastic emission mode for variable photomechanical effect to obtain pain information (analgesic effect)
4. Different applicators to adapt the therapy to the treatment area and deliver the correct dose in orthogonal mode
5. Pre-defined protocols



VIII. Spine pathologies

Treatment of low back pain with THEAL Therapy



VIII. Spine pathologies



18 ans, F, good health
Competitive gymnast,
9 training sessions per week
Low back pain with hyposthenia and dysesthesia of the left posteromedial thigh
L4-L5 disc herniation with 50% residual vertebral canal.

VIII. Spine pathologies

2 treatments per week for 5 weeks with Theal Therapy



6 months later

Reduction of the effect of compression and stenosis of the spinal canal
Stability of the clinical picture of a fully active athlete

Treatment examples – regenerative medicine

PRP



WOUNDS

I. PRP

- PRP is a type of regenerative medicine.
- New tissue formation and healing
- PRP combined with PBM accelerates the healing process.
- PBM with red laser light is the best combination with PRP.

Strengths of Mectronic Medicae lasers

1. Best wavelength range, emission modes and power to maximise PBM
2. High power range at 650 nm, from 10 mW to 2000 mW
3. Different applicators to adapt the therapy to the treatment area and to transfer the correct dose in orthogonal mode (Collimated or Aurix)




I. PRP

- PBM activates the microcirculation
- PBM with the 650nm wavelength acts directly on collagen
- PBM increases ATP production and accelerates the healing process
- PBM with high-powered red laser light is best for treating wounds

Strengths of Mectronic Medical lasers

1. The best range of wavelengths, emission modes and power
2. 2W 650nm with Theal and Chelt and up to 4W with iLux Plus
3. Temperature control system for better treatment efficiency and to activate the PBM in the best possible way (THEAL and CHELT)
4. Different applicators to adapt the therapy to the treatment area and to transfer the correct dose in orthogonal mode (e.g. ScanX system).
5. Predefined protocol



II. Wounds

Treatment of wounds with Theal Therapy



II. Wounds

Wound care: results after 15 treatments with Theal Therapy



Treatment examples – Podiatry

WARTS

ONYCHOMICOSIS

FOOT PATHOLOGIES

NAIL REGENERATION

I. Onychomycosis



BEFORE TREATMENT





1 MONTH AFTER






2 MONTHS AFTER

II. Warts



III. Nail regeneration



IV. Nail regeneration



Treatment examples – Wellness


SKIN REJUVENATION

I. Skin rejuvenation

- Anti-ageing, wrinkles and firmness
- PBM activates microcirculation
- PBM increases the production of collagen and elastin by photobiomodulation of the fibroblasts.
- The red laser light acts directly on the collagen.
- Non-invasive treatment. The epidermal layer is not damaged. This reduces the side effects of the treatment.

Strengths of Mectronic Medicales lasers

1. The best range of wavelengths, emission modes and power to maximise PBM.
2. 2W from 650nm with Theal and Chelt and up to 4W with iLux Plus
3. Temperature control system for improved treatment compliance and to activate PBM (Theal and Chelt only)
4. Different applicators to adapt the therapy to the treatment area and to transfer the correct dose in orthogonal mode.



I. Skin rejuvenation

Treatment of wrinkles with Theal therapy



Treatment examples - otorhinolaryngology

TINNITUS

I. Tinnitus

Treatment of tinnitus with Theal Therapy



Contraindications



The laser sources used emit in the visible and near infrared (NIR) spectrum.

There are two contraindications:

- Eye damage
- Thermal damage to the tissues exposed to the therapy, paying particular attention to the continuous emission and the phototype of the patient.

The first case is solved by the use of glasses, while the second is solved by the operator's ability to assess the thermal sensitivity and the correct phototype of the subject.

Avoid use in :

- **Patients with suspected or confirmed cancer**
 - **Pregnant women**
 - **Photosensitive areas**
 - **Infected areas**
 - **Epilepsy**
 - **Areas of bleeding**
- **Treatment in the sympathetic, vagal or cardiac areas in cardiac patients.**

PART FIVE: MECTRONIC LASER THERAPY DEVICES

I. Laser line 2023

MHEL



THEAL



CHELT



Laser line 2023

SMART

iLux SMART

Laser line 2023

PLUS

iLux PLUS

SMART

iLux SMART

Laser line 2023

XP

iLux XP

Ixyon XP

PLUS

iLux PLUS

SMART

iLux SMART

Laser line 2023

TOP

Ixyon XP

CHELT

XP

iLux XP

Ixyon XP

PLUS

iLux PLUS

SMART

iLux SMART

II. MHEL Therapy

MHEL THERAPY

Available models

iLUX SMART



iLUX PLUS



Available models

iLUX SMART



Available models

iLUX SMART



iLux SMART 810nm

7 W

iLux SMART 980nm

10 W

iLux SMART 1064nm

10 W • 15 W

iLux SMART 650nm

500 mW • 1 W

Available models

iLUX PLUS



Available models

iLUX PLUS



iLux PLUS 810+980nm

10 W • 15 W • 20 W • 30 W

iLux PLUS 1064nm

10 W • 15 W • 20 W • 30 W

iLux PLUS 810+1064nm

10 W • 15 W • 20 W • 30 W

iLux PLUS 650nm

2 W • 3 W • 4 W

Available models

	ILUX SMART	ILUX PLUS
AVERAGE POWER	7W – 15W	10W – 30W
PEAK POWER	Peak power : 20% more than the average power (for iLux 10W the peak power is 12W)	20W - 60W
NUMBER OF WAVELENGTHS	1	1 or 2
WAVELENGTH TYPE	<ul style="list-style-type: none">• 650nm• 810nm• 980nm• 1064nm	<ul style="list-style-type: none">• 650nm• 1064nm• 810nm + 980nm• 810nm + 1064nm

Available models

	ILUX SMART	ILUX PLUS
EMISSION MODE	10 EMISSION MODES <ul style="list-style-type: none">• CW• E²C• AntInf• PULS1• PULS2• PULS3• SINGLE PULSE• BURST• DIMMER• HPM	13 EMISSION MODES <ul style="list-style-type: none">• CW• E²C• AntInf• PULS1• PULS2• PULS3• SINGLE PULSE• BURST• DIMMER• HPM• THP• PBM• CUSTOM MODE
CUSTOM MODE	-	(WAVE EDITOR)

Available models

	ILUX SMART	ILUX PLUS
EFFECTS MODE	✓	✓
PATHOLOGICAL LIBRARY	60 protocols	150 protocols
CUSTOMIZED PROTOCOLS	✗	✓
INCLUDED APPLICATORS	1	2
OPTIONAL APPLICATORS	2	4
SCANX SYSTEM	✓	✓
SCREEN DIMENSIONS	7"	10,1"
FINGER SWITCH	✗	✓
FEEDBACK LED	✗	✓ during treatment

Available models

	ILUX SMART	ILUX PLUS
DIMENSIONS	295x265x147 mm	Without trolley: 295x265x150 Mm With trolley: 516x516x1017 Mm
TROLLEY	Not included	Included
WEIGHT	Without trolley : 3,5 Kg	Without trolley: 5 Kg With trolley: 24 Kg

The strengths of MHEL Therapy

Compact, lightweight
and portable

Affordable price for
a device with
excellent
performance

ScanX mode



Up to 30W power

More than 150
protocols

Handpiece with
finger switch

13 emission modes

III. THEAL Therapy

THEAL THERAPY

Available models

ILUX XP



IXYON XP



Available models

ILUX XP



Available models

ILUX XP



iLux XP 16 W

650 nm - 2 W



810 nm - 7 W



1064 nm - 7 W



iLux XP 30 W

650 nm - 2 W



810 nm - 8 W



980 nm - 10 W



1064 nm - 10 W



Available models

IXYON XP



Available models

IXYON XP



Ixyon XP 50 W	Ixyon XP 52 W
650 nm - 2 W	450 nm - 2 W
780 nm - 2 W	650 nm - 2 W
810 nm - 8 W	780 nm - 2 W
905 nm - 10 W	810 nm - 8 W
980 nm - 10 W	905 nm - 10 W
1064 nm - 10 W	980 nm - 10 W
	1064 nm - 10 W

COMING SOON

Available models

	ILUX XP 16W	ILUX XP 30W	IXYON XP 50W
AVERAGE POWER	16W	30W	50W
PEAK POWER	30W	50W	100W
MAXIMUM FLUENCE	30000 mJ/cm2	30000 mJ/cm2	50000 mJ/cm2
PULSE FREQUENCY	Up to 1000 Hz	Up to 1000 Hz	Up to 1000 Hz
PULSE ENERGY	Up to 30000 mJ	Up to 30000 mJ	Up to 50000 mJ
TIME GAP BETWEEN PULSES	1 – 1000 ms	1 – 1000 ms	1 – 1000 ms

Available models

	ILUX XP 16W	ILUX XP 30W	IXYON XP 50W
WAVELENGTHS NUMBER	3	4	6
WAVELENGTHS	650 + 810 + 1064 nm	650 + 810 + 980 + 1064 nm	650 + 780 + 810 + 905 + 980 + 1064 nm
MIX OF WAVELENGTHS	✓	✓	✓
LAMBDA WAVE CREATOR	X	X	✓
CALIBRATION SYSTEM	✓	✓	✓
THERMAL CONTROL	✓	✓	✓
VOBULATED THERMAL CONTROL	✓	✓	✓

Available models

	ILUX XP 16W	ILUX XP 30W	IXYON XP 50W
FINGER SWITCH	✓	✓	✓
EMISSION MODES	10	10	17
FEEDBACK LED	✓	✓	✓
NUMBER OF APPLICATORS	12	12	12
ZOOM APPLICATOR WITH THERMAL CONTROL	✓	✓	✓

Available models

	ILUX XP 16W	ILUX XP 30W	IXYON XP 50W
SPHERIC APPLICATOR	✓	✓	✓
AUTOMATIC MODE	✓	✓	✓
THERMAL CONTROL AUTOMATIC MODE	✓	✓	✓
NUMBER OF PROTOCOLS	More than 150	More than 150	More than 150
CUSTOM PROTOCOLS	✓	✓	✓
SCREEN DIMENSIONS	8"	8"	10,1"
POWER SUPPLY	100÷240V 50÷60Hz	100÷240V 50÷60Hz	100÷240V 50÷60Hz

Available models

	ILUX XP 16W	ILUX XP 30W	IXYON XP 50W
ABSORPTION	160 VA	160 VA	160 VA
DIMENSIONS	30 x 25 x 102 cm	30 x 25 x 102 cm	30 x 25 x 102 cm
WEIGHT	5 kg without trolley, 20 kg with trolley	5 kg without trolley, 20 kg with trolley	5 kg without trolley, 20 kg with trolley

The strengths of THEAL Therapy

Power from 10mW to 50w

Up to 6 different wavelengths

Patented wavelength mixing system

Up to 17 emission modes

Wave Creator

12 different applicators



Extensive interactive pathology library
(over 150 protocols)

Patented vobulated thermal control

Time, Joule and Trigger Point

Automatic working mode with
thermal control (ScanX)

Calibration system

	iLux Smart	iLux Plus	iLux XP	Ixyon XP	CHELT
Wounds	★ ★ (810 nm models) ★ ★ ★ (650 nm models)	★ ★ (810nm models) ★ ★ ★ ★ (650 nm models)	★ ★ ★ ★ ★	★ ★ ★ ★ ★	★ ★ ★ ★ ★
Podiatry	★ ★	★ ★ ★ ★	★ ★ ★ ★ ★	★ ★ ★ ★ ★	★ ★ ★ ★ ★
Pain management	★ ★ ★	★ ★ ★ ★	★ ★ ★ ★ ★	★ ★ ★ ★ ★	★ ★ ★ ★ ★
Acute pain	★ ★ ★	★ ★ ★ ★	★ ★ ★ ★ ★	★ ★ ★ ★ ★	★ ★ ★ ★ ★
Chronic pain	★ ★ ★	★ ★ ★ ★	★ ★ ★ ★ ★	★ ★ ★ ★ ★	★ ★ ★ ★ ★
Nerve inflammation, radiculopathy, neuropathy	★	★ ★	★ ★ ★ ★ ★	★ ★ ★ ★ ★	★ ★ ★ ★ ★

	iLux Smart	iLux Plus	iLux XP	Ixyon XP	CHELT
Muscular injuries	★ ★ ★	★ ★ ★ ★	★ ★ ★ ★ ★	★ ★ ★ ★ ★	★ ★ ★ ★ ★
Tendon injuries	★ ★ ★	★ ★ ★ ★	★ ★ ★ ★ ★	★ ★ ★ ★ ★	★ ★ ★ ★ ★
Oedemas	★ ★	★ ★ ★	★ ★ ★ ★ ★	★ ★ ★ ★ ★	★ ★ ★ ★
Active Trigger Points	★ ★ ★	★ ★ ★ ★	★ ★ ★ ★ ★	★ ★ ★ ★ ★	★ ★ ★ ★ ★
Latent Trigger Points	★ ★ ★	★ ★ ★ ★	★ ★ ★ ★ ★	★ ★ ★ ★	★ ★ ★ ★ ★
Fascia	★	★ ★	★ ★ ★ ★	★ ★ ★ ★	★ ★ ★ ★
Contractures	★ ★	★ ★ ★	★ ★ ★ ★ ★	★ ★ ★ ★ ★	★ ★ ★ ★ ★
Active and passive rehabilitation	★	★ ★	★ ★ ★	★ ★ ★	★ ★ ★
Physio-aesthetics			★ ★ ★	★ ★ ★	★ ★ ★ ★ ★